

PeakWorks Chromatography Software for any GC or LC

Software/Hardware Specifications

ADC	19 bit
Software	32 bit
Data Input Channels	Detector Input -2; optional- 4 Temperature Input- 2
Timed events	Three; optional- 6
Communication	USB port; optional converter available
PC Required	Pentium, 32 Mbytes RAM, Windows 95, 98 or NT
Reintegration	Yes
Reports	Rich text, PRN, Printer
Calibration	Multipoint
Auto run	Valves can be setup to automatically inject samples over long periods of time